

採尿方法

How to take urine

別紙4

Separate Sheet No.4

1) 採尿カップに尿を取る

Take urine in a
test cup.



2) 採尿容器の腹を 押さえ、採尿カップの 尿を吸い上げる

Pushing the side of the
bottle,
pump up urine from
the cup by a bottle.



3) 線のところまで 吸い上げ、容器の ふたを閉める

Fill the bottle to a
marker line, and
screw the cap
back on.



4) 提出用袋に 所属と氏名を記入 し提出

Write your
affiliation and name
on the bag.

